

# Portable Wired Ultrasonic Plus LED Skin Care Machine 3 Colors



Spa & Equipment SKIN+ACT

Toll Free: 1.800.418.0830	<a href="http://www.skinact.com">www.skinact.com</a> <a href="http://www.spaandequipment.com">www.spaandequipment.com</a>	International: 818.834.6640
---------------------------	--	-----------------------------

## NOTICE

No part of this document may be photocopied, reproduced, or translated into any other language without prior written consent of Spa and Equipment/SkinAct.

The information in this document is subject to change without notice at any time and should not be construed as a commitment by Spa and Equipment/SkinAct.

Spa and Equipment/SkinAct assumes no responsibility for any errors that may be in this document, nor does it make expressed or implied warranty with regard to this material, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

Spa and Equipment/SkinAct shall not be liable for incidental or consequential damages arising out of the furnishing, performance, or use of this document and the program material which it describes.

Spa & Equipment SKIN+ACT




Toll Free: 1.800.418.0830

[www.skinact.com](http://www.skinact.com)  
[www.spaandequipment.com](http://www.spaandequipment.com)

International: 818.834.6640

# Measurements



Spa & Equipment SKIN+ACT

# Highlights

## Features:

3Mhz vibrations infiltrate the dermis and promote recovery of elasticity and collagen so that fine wrinkles disappear.

This powerful photon therapy uses the natural healing power of light to treat numerous skin issues.

Utilizes the latest in LED technology to ensure that you get the best treatment as possible.

Wired operated (not operated by battery)

## Benefits:

Brightening

Improve blood circulation

Anti-wrinkle

Balance skins and enhance blood circulation

Remedy acne and dermatitis

During a treatment, the user feels no heating, no strong vibrations but feels soft and small vibrations as the Micro Vibration waves frequency is high and the skin cannot feel much. But pigment gets broken down, molecules get activated, and blood circulation and metabolism gets motivated during treatment. Skin becomes finer and smoother, pore size shrinks, wrinkles released, and pigmentation fades away.

Blue LED light penetrates about 0.5mm. It helps calm the skin and fight skin rashes; suitable for the skin which is sensitive, oily, and with skin rashes. It can adjust skin sebum secretion. It also helps to get rid of acne. Blue light also helps with skin tightening.

Red LED light penetrating 8-10mm is suitable for any type of skin; it can reach the derma to stimulate fibre cells. Warming can improve blood circulation and collagen albumen hyperplasia to reduce wrinkles and make skin more elastic and smooth. The mild wavelength is suitable for around the eyes and forehead.

Green LED light between red light and blue light is suitable for the mixed skin. It can improve the function of cells and oxygen and activate the microcirculation function of body surface. It also helps to calm the skin and make it smoother.

Specifications: Power source: Input: AC100V-240V, 50/60Hz, 0.25A; Output: DC15V 0.4A; Machine net weight: 150g; Package cross weight: 380g; Dimensions: 21\*15\*6cm.