

Rover Medical Treatment Table Specifications Sheet



Spa & Equipment SKINACT

Toll Free: 1.800.418.0830

www.skinact.com
www.spaandequipment.com

International: 818.834.6640

NOTICE

No part of this document may be photocopied, reproduced, or translated into any other language without prior written consent of Spa and Equipment/SkinAct.

The information in this document is subject to change without notice at any time and should not be construed as a commitment by Spa and Equipment/SkinAct.

Spa and Equipment/SkinAct assumes no responsibility for any errors that may be in this document, nor does it make expressed or implied warranty with regard to this material, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

Spa and Equipment/SkinAct shall not be liable for incidental or consequential damages arising out of the furnishing, performance, or use of this document and the program material which it describes.





Toll Free: 1.800.418.0830

www.skinact.com
www.spaandequipment.com

International: 818.834.6640

Measurements




Measurements



KIN|CT

Spa & Equipment

Measurements



SKINACT

Spa & Equipment

Highlights

Movements:

- Electric height adjustment
- Manual backrest adjustment (air-gas pump assisted)
- Manual leg rest adjustment (air-gas pump assisted)
- Manual rolling wheels for easy table movement

Dimensions/Specifications:

- Length: 78"
- Width: 27"
- Heights: 22.5" – 36"
- BASE ONLY Dimensions: 59.5"L x 24"W
- Backrest angle: 0° flat to 79° incline
- Leg rest angle: 0° flat to 28° incline - 0° flat to 49° decline
- Weight: 245 lbs
- Capacity of 600 lbs
- Color available: Grey
- Ultra-durable super fiber PU upholstery with a lifespan of over 10 years. (It's the same quality upholstery used in automobile manufacturing)
- Easy to move between rooms with the lockable wheels
- Heavy-duty build makes the base super stable
- Table paper holder included
- Ultra-quiet motor

* All height measurements are from the ground to the top of the cushion *

** All measurements have a 0.25 to 0.75 inch tolerance **

Warranty:

- 1-year warranty

