# Bliss Spa Facial Treatment Table Specifications Sheet





Toll Free: 1.800.418.0830

www.skinact.com www.spaandequipment.com

International: 818.834.6640

#### NOTICE

No part of this document may be photocopied, reproduced, or translated into any other language without prior written consent of Spa and Equipment/SkinAct.

The information in this document is subject to change without notice at any time and should not be construed as a commitment by Spa and Equipment/SkinAct.

Spa and Equipment/SkinAct assumes no responsibility for any errors that may be in this document, nor does it make expressed or implied warranty with regard to this material, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

Spa and Equipment/SkinAct shall not be liable for incidental or consequential damages arising out of the furnishing, performance, or use of this document and the program material which it describes.














Toll Free: 1.800.418.0830

www.skinact.com www.spaandequipment.com

International: 818.834.6640

# **Measurements**






# **Measurements** 44" 55" 70° Backrest Angle 28" \KIN\CT Spa\text{Equipment}

# **Wheels and Controls**







# **Highlights**

#### **Movements:**

- Electric backrest angle adjustment
- Electric height adjustment
- Manual rolling wheels for easy table movement

### **Dimensions/Specifications:**

- Length: 73.5"
- Width: 28"
- Heights: 21" 34"
- Base dimensions (BASE ONLY): 55" x 24"
- Backrest angle: 0° flat to 70° incline
- Capacity of 500 lbs
- Cushion thickness: 4"
- Color available: White
- · High-quality PU upholstery; ultra-durable material
- Removable pillow included
- Breather hole for when clients are face down
- Four wheels
- Heavy-duty construction
- · Control with hand remote or with buttons on the bed itself
- \* All height measurements are from ground to top of cushion \*
- \*\* All measurements have a 0.25 to 0.75 inch tolerance \*\*

#### **Warranty Information:**

1-year warranty

